
Apto® F 01 Full Face Mask

Not for sale in the US. Pending 510(k) clearance.

English User Instructions	1
Dansk Brugervejledning	7
Deutsch Bedienungsanleitung	13
Ελληνικά Οδηγία χρήσης	19
Español Instrucciones de uso	25
Suomi Käyttöohjeet	31
Français Instructions d'utilisation	37
Italiano Istruzioni per l'utente	43
日本語 ユーザー取扱説明書	49
Bahasa Melayu Arahan Pengguna	55
Norsk bokmål Brukerveiledning	61
Nederlands Gebruikershandleiding	67
Polski Instrukcja obsługi	73
Português Instruções do utilizador	79
Srpski Korisničko uputstvo	85
Svenska Bruksanvisning	91
Türkçe Kullanım talimatı	97
العربية إرشادات المستخدم	103
עברית הוראות למשתמש	109

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1 Introduction

1.1 Intended Use

The Apto® F01 Full Face Mask is intended for adults weighing >30kg (66lbs) who have been prescribed PAP or Bi-Level therapy. The masks are intended for single patient use in the home and/or multi-patient use in the hospital or other clinical setting where proper disinfection of the device can occur between patient uses.

1.2 Symbols

	Product number 009972: Apto® F01 FFM Small 009973: Apto® F01 FFM Medium 009974: Apto® F01 FFM Large		
	Batch code		Keep away from sunlight
	Medical device		Single patient, multiple use
	Date of manufacture		Expiration date
	Manufacturer information		Temperature Limitation
	Read user instructions		Atmospheric pressure limitation
	Swiss authorized representative		Humidity Limitation
Rx Only	The Apto® F01 Full Face Mask shall be prescribed by a medical professional (where applicable) or in accordance with local regulations. (Applies to the U.S. only)		
	Apto® F01 Full Face Mask meets all requirements for CE marking according to relevant European health, safety and environmental protection legislation.		
	Apto® F01 Full Face Mask meets all requirements for UKCA marking according to applicable United Kingdom health, safety and environmental protection legislation.		

1.3 Corporate Information

Manufacturer	Breas Medical AB Företagsvägen 1, SE-435 33 Mölnlycke, SWEDEN
UK REP	Breas Medical Ltd. Unit A2, Timothy's Bridge Road, Stratford-upon-Avon, Warwickshire, UK

Apto, Breas Apto, Breas and the Breas logo are trademarks of Breas Medical AB.

1.4 Mask Parts (Figure A)

See *Figure A: Mask Parts test*, page 115.

1. Elbow, with integrated Anti-asphyxia valve (AAV)	4. Headgear clips
2. Frame	5. Headgear
3. Cushion, with integrated leak port	

2 Safety Information

The mask shall only be used with compatible PAP or Bi-level equipment recommended by your clinician, and with settings recommended by the clinician. Using the Apto® F01 Full Face Mask in combination with incompatible medical devices can have adverse effects on the safety and performance of the mask.

- The mask may not be suitable for those predisposed to aspiration.
- Consult a physician or dentist if you encounter tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.
- This mask is not suitable for providing life support ventilation.
- If a patient has any of the following conditions, the prescribing clinician shall decide if the use of a full-face mask is the best option or if a different interface is preferred: Facial burns, recent head or facial surgery.
- Do not use the mask with an open skin wound. Consult your prescribing clinician for advice.
- Mask use with oxygen: When your PAP or Bi-level device is not operating, the oxygen flow must be turned off. Otherwise, oxygen might accumulate in the device enclosure and create a risk of fire. You should not smoke or be in the vicinity of someone smoking while using oxygen with your mask.

WARNING!

Apto® F01 Full Face Mask is a full face mask with a leak port for minimising re-breathing and an anti-asphyxia valve (AAV) for allowing spontaneous breathing in case of failure of the air flow to the mask.

- The mask shall not be used without a continuous flow of air.
- Do not block the leak port or the anti-asphyxia valve (AAV). Occlusion of the leak port leads to rebreathing of exhaled gases and adverse effect on the safety and quality of the therapy.
- For patients unable to remove the mask themselves, the mask must only be used under qualified supervision and by advice from a clinician.

If you experience any adverse reaction to any part of the mask system, stop using it and consult your physician or provider. Adverse reactions may include, but are not limited to:

- Discomfort or irritation in or near the eyes or anywhere on the head.
- Excessive leaks between mask seal and face.
- Gum, jaw or mouth soreness/ discomfort or an aggravation of an existing dental condition.
- Allergic reaction, mask discomfort or irritation.
- Nausea or vomiting.

NOTE

Any serious incident that has occurred in relation to this device should be reported to the competent authority and the manufacturer.

3 Operating Instructions

3.1 Checks Before Use (Figure B)

See *Figure B: Checks Before Use*, page 116.

- 1 Inspect for damage or deterioration. If there is visible damage (cracks, tears, etc.), do not use your mask until the damaged part is replaced.
Check that the mask is clean and free from objects that may prevent breathing.
- 2 Check that all parts of the mask are assembled, see *Assembling the Mask (Figure F)*, page 5.
- 3 Attach one end of the patient circuit to the elbow and the other end to your PAP or Bi-level device.
- 4 Start the treatment on the PAP or Bi-level device and hold the mask to your face. Check that air flows through the leak port. If there is no flow, do not use the mask.
- 5 Stop the treatment and hold the mask to your face. Check that you can breathe through the Anti-asphyxia valve. If you cannot breathe through the valve, do not use the mask.

3.2 Putting on Your Mask (Figure C)

See *Figure C: Putting on the Mask*, page 117.

Before putting on your mask, perform the *Checks Before Use (Figure B)*. The patient circuit may be connected to the mask either before or after putting the mask on.

- 1 Put the mask on your face with the cushion covering your mouth and nose. Pull the headgear over your head.
- 2 Attach the clips to the frame.
- 3 The mask should have a leak-free yet comfortable fit. If required, adjust the fit using the strap tabs. First, adjust with the upper tabs and then with the lower tabs. If a good fit cannot be accomplished, try a different size of headgear or cushion, see section 6 *Spare Parts*.

To remove the mask: first, rotate the clips to release them from the frame, and then pull the headgear over your head.

4 Cleaning

WARNING!

Using the mask with multiple patients without proper disinfection may cause cross-contamination.

The cleaning procedure described in this section applies to single patient use. For multiple patient use, the mask shall be disinfected by a healthcare professional between patients, see 7.2 *Reprocessing*.

WARNING!

Risk of reduced performance.

Clean as described in this manual only. The use of other intervals, methods, or agents may compromise the function of the mask and consequently the safety or the quality of the therapy.

Cleaning the Parts (Figure D)

See *Figure D: Cleaning*, page 118.

This cleaning procedure applies to single patient use.

- 1 Disassemble the parts to clean.
Daily cleaning: Cushion and elbow
Weekly cleaning: Frame and headgear.
- 2 Clean the parts with warm water (20 to 40°C, 68 to 104°F) and a mild liquid detergent, using a soft bristle brush (like a toothbrush). Do not use disinfectants, antibacterial soaps, scented solutions, moisturisers or oils to clean the parts. Check that the parts are visually clean before proceeding.
- 3 Rinse the parts thoroughly with water, for at least 15 seconds. Check that all detergent residues are rinsed before proceeding.
- 4 Allow the parts to air dry, out of direct sunlight, for at least 30 minutes. Check that the parts are completely dry before proceeding.
- 5 Inspect that the parts are undamaged, visually clean and dry before assembling the mask. If required, repeat the relevant cleaning steps. If a part is damaged or cannot be cleaned completely, replace the part before using the mask.

Disassembling the Mask (Figure E)

See *Figure E: Disassembling the Mask*, page 119.

- 1 Detach the headgear clips by rotating them.
- 2 Slide the upper headgear loops from the frame hooks.
- 3 Detach the elbow from the frame by squeezing the elbow latches and pulling the elbow out.
- 4 Detach the cushion from the frame by pulling at the rigid part that connects to the frame. Avoid pulling on the soft part of the cushion.

Assembling the Mask (Figure F)

See *Figure F: Assembling the Mask*, page 120.

- 1 Place the cushion in your palm with the frame on top of it. Tap the frame to attach the parts (indicated by a clicking sound).
- 2 Push the elbow onto the frame so it attaches (indicated by a clicking sound).
- 3 Slide the upper headgear loops into the frame hooks, ensuring that the Breas logo is on the outside of the headgear.

5 Disposal

Dispose of the mask according to local environmental regulations. If required, first disassemble the parts.

6 Spare Parts

Part	Breas Article Number
Cushion	Small: 009975 Medium: 009976 Large: 009977
Mask Frame	009982
Headgear	Small: 009978 Medium: 009979 Large: 009980
Headgear Clips (2 pieces)	009981
Mask system without headgear (US only)	Small: 009983 Medium: 009984 Large: 009985
Elbow	009986

7 Technical Information

Specification	Value
Patient circuit connector	22 mm male (ISO 5356-1)
Operating conditions	Temperature: 5–40 °C (41–104 °F) Humidity: 15 – 95 %RH
Expected lifetime (single patient, multiple use)	1 year
Storage and transport conditions	-20 – 60 °C (-4 – 140 °F), clean and dry, protected from sunlight.
Shelf life	2 Years
Pressure range For pressure-flow diagram, see <i>Figure G</i> .	4–40 cmH ₂ O
Resistance	@ 50 l/min: < 0.2 cmH ₂ O @ 100 l/min: < 0.6 cmH ₂ O
Deadspace	L: 311 ml, M:281 ml, S:243 ml
AAV system pressure for opening/closing to atmosphere	< 2.4 cmH ₂ O
Resistance with AAV open to atmosphere (@50 l/min)	Inspiration: < 0.6 cmH ₂ O Expiration: < 0.7 cmH ₂ O
Sound power level	≤ 26 dB(A), with 1 dB(A) measurement uncertainty.
Sound pressure level @1m	≤ 20 dB(A), with 1 dB(A) measurement uncertainty.

7.1 Materials

The Apto® F01 Full Face Mask is not made with natural rubber latex. It is not made with hazardous materials and is compliant with applicable parts of ISO 10993:1, ISO 18562:1 and the REACH regulation (EU).

7.2 Reprocessing

The Apto® F01 Full Face Mask can be reprocessed for multiple patient use. The reprocessing is to be performed by qualified healthcare professionals only.

For reprocessing instructions, see <https://www.breas.com/apto>.

8 Frequently Asked Questions

For additional product information, see <https://www.breas.com/apto>.

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