

Breas Apto® N 01

Technical data

Breas Apto® N01													
Pressure/Flow Curve													
Pressure flow diagram <table border="1"><caption>Pressure flow diagram data points (approximate)</caption><thead><tr><th>Pressure (cmH₂O)</th><th>Flow rate (l/min)</th></tr></thead><tbody><tr><td>4</td><td>18</td></tr><tr><td>11</td><td>30</td></tr><tr><td>17</td><td>38</td></tr><tr><td>24</td><td>45</td></tr><tr><td>30</td><td>52</td></tr></tbody></table>		Pressure (cmH ₂ O)	Flow rate (l/min)	4	18	11	30	17	38	24	45	30	52
Pressure (cmH ₂ O)	Flow rate (l/min)												
4	18												
11	30												
17	38												
24	45												
30	52												
Therapy Pressure Range	4-30 cm/H ₂ O												
Dead Space	L: 96 ml, M: 91 ml, S: 81 ml												
Resistance	@50 l/min: < 0.3 cm/H ₂ O @100 l/min < 0.13 cm/H ₂ O												
Sound Pressure Level @ 1m	< 17 dB(A)												
Dimension (packaging) H x W x D	155 x 125 x 110 mm												
Weight (mask)	<82 g												
Weight (Inc packaging)	<250 g												
Shelf Life	5 Years												
Expected Service Life	1 Year												
Warranty	90 Days												
Frame Material	Nylon												
Headgear Material	Nylon/Elastane												
Headgear Clips Material	Nylon												
Cushion Material	Silicon/Polycarbonate												
Elbow Body	Polycarbonate												
Swivel	Nylon												
Elbow Latch	Nylon												
Storage and Transport Conditions	-20– 60 °C (-4– 140 °F), clean and dry protected from sunlight												
Operating Conditions	Temperature: 5–40 °C (41–104 °F) Humidity: 15– 95 %RH												

Be more™

References: Breas Apto® N01 Instructions for Use; Breas Data on File.
Breas, the Breas logo, Apto, and Breas Apto are trademarks of Breas Medical AB. © 2026 Breas Medical AB – All rights reserved. Breas Medical AB reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information.
Issued: May 2026. MAR 008208 Rev_0



Breas Medical AB
Företagsvägen 1, SE-435 33
Mölnlycke, Sweden
+46 31 86 88 00



For more
information
please visit
breas.com/apto