



# Vivo 45LS

by Breas

## Prescription Form

### Patient Information

|               |             |                          |                                                                     |
|---------------|-------------|--------------------------|---------------------------------------------------------------------|
| Patient Name: | Ventilator: | Vivo 45LS                | Do not substitute                                                   |
| Phone:        | DOB:        | Estimated length of need | months (99=lifetime)                                                |
| Diagnosis:    | Date:       | Hours of use:            | During Sleep PRN                                                    |
|               |             | HCPCS Code:              | E0465 Home ventilator, invasive E0466 Home ventilator, non-invasive |

### Profile 1 Mode

|                         |                                                                       |                                          |                                          |                       |                        |          |
|-------------------------|-----------------------------------------------------------------------|------------------------------------------|------------------------------------------|-----------------------|------------------------|----------|
| Modes:                  | PSV                                                                   | PSV(TgV)                                 | PSV(TgV)+AE                              | PCV(A)                | PCV(A+TgV)             | PCV-SIMV |
|                         | VCV(A)                                                                | VCV-SIMV                                 | CPAP                                     | PCV-MPV               | VCV-MPV                | HFNT     |
| Insp Pressure (cmH2O)   | Pressure Support (cmH2O)                                              | PEEP (cmH2O)                             | CPAP (cmH2O)                             |                       |                        |          |
| Resp Rate (bpm):        | Insp. Time (s):                                                       |                                          |                                          |                       |                        |          |
| Tidal Volume (mL):      | HFNT (lpm):                                                           |                                          |                                          |                       |                        |          |
| Target Volume:          | Add Target Volume (with applicable mode)                              | Target Volume (mL):                      | Max Pressure (cmH2O):                    | Min Pressure (cmH2O): |                        |          |
|                         | Insp. Pressure                                                        |                                          |                                          |                       |                        |          |
| Auto-EPAP:              | Add Auto-EPAP                                                         | PC Breath                                |                                          |                       |                        |          |
| Auto-EPAP Settings:     | Max PS (cmH2O)                                                        | Min PS (cmH2O)                           | Max EPAP (cmH2O)                         | Min EPAP (cmH2O)      | Pressure Limit (cmH2O) |          |
| Sigh Breath:            | Add Sigh Breath (with applicable mode)                                | Sigh Rate:                               | Sigh Rate: Sigh % of Volume or Pressure: |                       |                        |          |
| Ramp on/off (10-60 min) | Starting Pressure                                                     |                                          |                                          |                       |                        |          |
| O2 (LPM):               |                                                                       |                                          |                                          |                       |                        |          |
|                         | Titrate Rise Time, Insp. Trigger and Exp. Trigger for patient comfort | +/- 100 mL volume adjustment (if needed) | Add Home Adjust Settings                 |                       |                        |          |

Notes:

### Monitoring

|                                         |            |                                       |         |            |   |     |
|-----------------------------------------|------------|---------------------------------------|---------|------------|---|-----|
| Monitor via Vivo integrated monitoring: | SpO2       | EtCO2                                 | PtcCO2  | Continuous | Q | Hrs |
| Download Data Every                     | Days       | Enroll in EveryWare remote monitoring |         |            |   |     |
| Apnea Interval:                         | 15 seconds | Other                                 | seconds |            |   |     |

### Additional Profiles

|                   |               |          |          |              |               |
|-------------------|---------------|----------|----------|--------------|---------------|
| Profile 2 - Mode: | Press(cmH2O): | Vol(mL): | RR(bpm): | PEEP(cmH2O): | O2 Flow(lpm): |
| Notes:            |               |          |          |              |               |
| Profile 3 - Mode: | Press(cmH2O): | Vol(mL): | RR(bpm): | PEEP(cmH2O): | O2 Flow(lpm): |
| Notes:            |               |          |          |              |               |

### Physician Information

|                      |       |
|----------------------|-------|
| Physician Name:      | NPI#: |
| Physician Signature: | Date: |

Vivo by Breas, Breas, EveryWare by Breas, eSync and the Breas logo are trademarks of Breas Medical AB. © 2026 Breas Medical AB - All rights reserved. Breas Medical reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information.

MAR-008315 Rev 1



Breas Medical, Inc.  
16 Esquire Rd.  
N. Billerica, MA 01862  
+1 855 436 8724



For more  
information  
visit  
**breas.us**